




Clinical  
Supplies  
Management  
Team



## As of now the new blue topped transfusion sample tubes must be used for all paediatric transfusion samples

New blue topped paediatric blood transfusion sample tubes are now available to order via NDC. All ward top up/order books will be updated accordingly.

For babies and infants **up to 10kg** please use the blue top 1.2ml tube - order sku 185641 product code 06.1662.002



For infants and children **over 10kg** please continue to use blue top 4.5ml tube - order sku 183005 product code 05.1554.009 (this is the same as the adult transfusion sample tube)

**Please note:** Pink topped Teklab 1.3ml TLK 200PP will **no longer be acceptable for transfusion samples** (but must continue to be used for haematology requests please)

Any queries please contact Carol Balfour, Procurement 0131 465 7656 or Jane Oldham, Transfusion Practitioner SJH 53164